

CPR stat•padz® & CPR-D•padz®

Benefits of Real CPR Help™ Come with
Placing Electrodes



Get Real CPR Help™ on ZOLL Defibrillators
by Placing Electrodes. There are No
Extra Sensors, Steps or Cables Required.

ZOLL
Advancing Resuscitation. Today.™

Real CPR Help

Failure to adequately compress a victim's chest is a common error during CPR.^{1,2,3} Real CPR Help converts compression data instantaneously, and provides visual and audio prompts to help the rescuer with the proper rate and depth. Activated audio prompts say, "Push Harder" or "Good Compressions."

See-Thru CPR

See-Thru CPR^{4,5} is an advanced feature on ZOLL's R Series[®], E Series[®], and AED Pro[®]. Artifact caused by CPR compressions is minimized giving clinicians the ability to see organized electrical activity during compressions without interruptions.

CPR Index

The CPR Index⁶ is a feature on ZOLL's R Series and E Series. It is an intuitive visual summary of chest compression efforts that takes into account rate, depth and gaps between CPR cycles.

CPR-D•padz

The **CPR-D•padz** was designed for ZOLL AEDs. The unique one-piece design uses a CPR landmark for hand placement, and makes it easy to place the **CPR-D•padz** on the victim quickly. The **CPR-D•padz** also lasts five years – the longest life of any electrodes sold today.



CPR stat•padz

CPR stat•padz is designed for use with ZOLL manual defibrillators. The two piece construction ensures all therapy options are available to the ALS responder. The unique design gives it utility in both the Anterior/Anterior and Anterior/Posterior placement.



	Shelf Life	Construction	Defibrillate	Pace	Pre-Connect	Real CPR Help	See-Thru CPR ^{4,5}	CPR Index ⁶
CPR-D•padz	5 Yrs	One Piece	■		■	■	■	■
CPR stat•padz	2 Yrs	Two Piece	■	■	■	■	■	■

General Specifications

	CPR stat•padz	CPR-D•padz
Conductive Gel:	Solid, adhesive polymer	Solid, adhesive polymer
Conductive Element:	Tin	Tin
Packaging:	Multilayer laminate	Multilayer laminate
Length, Lead Wires:	49 in (124.5 cm)	36 in (91cm)
Apex Electrode		
Length:	6.5 in (16.5 cm)	6.1 in (15.5 cm)
Width:	4.0 in (10.2 cm)	5.6 in (14.1 cm)
Total Area:	26 in ² (168.3cm ²)	34.1 in ² (86.6 cm ²)
Length, Conductive Gel:	5 in (12.7 cm ²)	3.5 in (8.9 cm)
Width, Conductive Gel:	3.5 in (8.9cm)	3.5 in (8.9 cm)
Area, Conductive Gel:	17.5 in ² (113 cm ²)	12.3 in ² (79 cm ²)
Sternum Electrode		
Length:	6.5 in (16.5 cm)	6.1 in (15.5 cm)
Width:	4.0 in (10.2 cm)	5 in (12.7cm)
Total Area:	26 in ² (168.3cm ²)	30.5 in ² (77.5cm ²)
Length, Conductive Gel:	5 in (12.7 cm)	3.5 in (8.9 cm)
Width, Conductive Gel:	3.5 in (8.9cm)	3.5 in (8.9 cm)
Area, Conductive Gel:	17.5 in ² (113 cm ²)	12.3 in ² (79 cm ²)
Central Section		
Length:	N/A	12.9 in (32.9 cm)
Width:	N/A	4.2 in (10.6 cm)
Length, CPR Device:	3.4 in (8.7cm)	3.4 in (8.7cm)
Width, CPR Device:	2 in (5.1 cm)	2 in (5.1 cm)

¹ Moser DK, Dracup K, Guzy PM, Taylor SE, Breu C. Cardiopulmonary resuscitation skills retention in family members of cardiac patients. American Journal of Emergency Medicine. 1990;498-503.

² Kern KB, Hilwig RW, Berg RA, Ewy GA. Efficacy of chest compression-only BLS CPR in the presence of an occluded airway. Resuscitation.1998;39:179-188.

³ Handley AJ, Handley JA. The relationship between rate of chest compression and compression:relaxation ratio. Resuscitation. 1995;30:237-241.

⁴ 510(k) pending

⁵ Defibrillator dependent.

Currently See-Thru CPR is pending 510K clearance for the AED Pro and pending Health Canada medical Device License for the AED Pro and E Series.

Specifications subject to change without notice.

© ZOLL Medical Corporation 2007. All rights reserved.

ZOLL, **CPR-D•padz** and **CPR stat•padz** are registered trademarks of ZOLL Medical Corporation. Advancing Resuscitation. Today., R Series, Real CPR Help, CPR Index and See-Thru CPR are trademarks of ZOLL Medical Corporation.

Printed in U.S.A. 050710 9656-0199

ZOLL Medical Corporation Worldwide Headquarters

269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-348-9011

ZOLL Direct Sales Subsidiaries

Australia, +61 2 8424 8700
www.zoll.com.au

Austria, +43 650 413 6222
www.zollmedical.at

Canada, 905-629-5005
www.zoll.com

France, +33 (1) 30 05 14 98
www.zoll.fr

Germany, +49 2236 87 87 0
www.zollmedical.de

The Netherlands, +31 (0) 488 41 11 83
www.zoll.nl

United Kingdom, +44 (0) 1928 595 160
www.zoll.com

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.

