

Swaddler for Phototherapy

Bionix® Swaddler for Phototherapy gives care givers the ability to safely swaddle infants undergoing phototherapy, while allowing approximately 85% light transmittance.



PROBLEM

During the first week of life, approximately 1 in 8 full term infants are readmitted for jaundice, and nearly all premature infants weighing <1500g become clinically jaundice. Add bright lights, sticky leads, and various other attachments, and the infant is suddenly in a high stress environment.

SOLUTION

Swaddler for Phototherapy provides a familiar barrier for the arms and legs that does not interfere with phototherapy treatment or infant monitoring. The unique fabric allows for approximately 85% of light transmittance to reach the infants skin.

RESULTS

A simple solution that provides a comforting barrier to infants undergoing treatment of hyperbilirubinemia.

PRODUCT OFFERING

#1100 - 10 Premie (<1500g) Swaddlers

#1105 - 10 Newborn Regular (>1500g) Swaddlers

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

Photo Therapy Light Source	Average Light Transmittance
GE Healthcare Giraffe® SPOT PT Lite™ Phototherapy System	86.10%
Dräger Photo-Therapy 4000	83.06%

Note: An Olympic Model 22 Bili-Meter™ with type B-22 sensor was used for the light transmission testing. Giraffe, and Spot PT Lite are trademarks of Datex-Ohmeda, Inc.



EC REP Obelis s.a • Boulevard Général Wahis 53 • 1030 Brussels, BELGIUM
 Tel: +(32) 2.732.59.54 • Fax: +(32) 2.732.60.03 • E-Mail: mail@obelis.net
 Bionix LLC • 1670 Indian Wood Circle, Maumee, Ohio 43537, United States
 Phone: +1 419.727.8421 • Fax: +1 800.455.5678 • www.Bionix.com



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • www.bionix.com

©2021 Bionix LLC | US Patent # 8,756,731 | MKT-0048 | Rev. C