

AUGUSTINE

TREATMENT RECLINER

PREMIUM COLLECTION



Developed in conjunction with Oncology clinics, the Augustine was designed to be 100% patient-centric; to be comfortable and help put the patient at ease.

- Enhanced range of motion for optimum, patient selectable, positioning
- Spacious seat provides personal space and luxurious comfort
- Thoughtful features to ease chair movement and cleaning
- Unmatched stability and refined styling with concealed casters

Augustine demonstrates Winco's commitment to developing products that meet the needs of the people that live in and touch the recliners on a daily basis.

**ENHANCING THE PATIENT
AND CAREGIVER EXPERIENCE.**



SERENITY: PERSONAL SPACE & INDEPENDENCE

Often the medical treatment environment unavoidably impinges upon the patient's personal space. Augustine addresses this with:

- A low seat height to ease entry and exit
- A spacious 26" seat width for extra room to stretch out
- Easy access recline controls that allow patients to obtain desired positioning without assistance

When asked about their needs, patients were clear they wanted a chair that provided them with independence and a relaxing experience.

AUGUSTINE IS DESIGNED TO MIMIC FURNITURE THAT YOU WOULD BE PROUD TO HAVE IN YOUR OWN HOME



The Central Front Braking option

CLINICIAN INPUT: DESIGN ENHANCEMENT

Facility managers requested an easier way to thoroughly clean the recliner. Once removed, the Quick-Release Seat provides full access to the inner sides of the recliner, backrest and seat for fast and thorough cleaning.

Simply pull-up at the front and slide forward. The seat is securely tethered so it won't get separated from the recliner. To replace simply line up the tabs, slide into place and push down – done!

QUICK RELEASE SEAT: NO NEED FOR ANY TOOLS AND RISK OF LOST PARTS



CLINICAL UTILITY: PRACTICAL INNOVATION

Augustine has the clinical features required in a high use medical environment – casters for ease of movement, rigorously tested for durability and thoughtful features for ease of cleaning. These chairs bring dignity, comfort and grace back to the clinical setting without forfeiting the necessary clinical utility.

- Standard rear central caster lock for impressive stability
- Quick Release Seat is easily removed for efficient and thorough cleaning
- An optional Central Front Brake allows secure braking from the front of the chair, with no obtrusive levers or mechanisms to catch on or impede staff, patients or visitors.



ANATOMY OF THE AUGUSTINE TREATMENT RECLINER



SPECIFICATIONS

MODEL	Augustine Treatment Recliner (5001)
Arm Style	Fixed
Casters	Front: 3" Fixed Rear: 3" Swivel Central Locking
Weight Capacity	500 lbs (226.80 kg)
Overall Height	45.75" (116.2 cm)
Overall Width*	34.5" (87.63 cm)
Seat Height (from Floor)	19.75" (50.2 cm)
Seat Width (Between Arms)	26" (66.04 cm)
Seat Depth	20.75" (52.71 cm)
Seat Back Height	28.5" (72.39 cm)
Armrest Height	26.75" (67.9 cm)
Upright Length	40.0" (101.6 cm)
Recline Length	76.0" (193.00 cm)
Range of Seat Angle	75° (Upright) 12° (Reclined)
Weight	181 lbs (82.1 kg)

OVER 50
VINYL OPTIONS
AVAILABLE



* Add 2.5" to overall width for each table in the down position.

* Add 11" to overall width for each raised side table.

All specifications subject to change without notice, please contact our customer care team to verify any critical dimensions prior to ordering.

SUGGESTED OPTIONS & ACCESSORIES



IV Pole (IV)

Left, Right or Both Sides



Accessory Hook (AH)

Left, Right or Both Sides



Heat Only (HT)



Heat & Massage (HM)



Premium Side Table (PT)

Left, Right or Both Sides



Pivot Table (PV)

Left, Right or Both Sides



Central Front Brake (BE)

WARRANTY Contact Winco or visit our website for warranty details. **COLORS** Contact Winco or visit our website for all available color selections. **NOTE** Winco reserves the right to change specifications without notice. **CERTIFICATIONS** 1) The Augustine Treatment Recliner meets ANSI/BIFMA X5.4 for reliability testing. 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture. 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for upholstered furniture. 4) Heat and Massage and Heat Only Options meet electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92.



CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com