

# Lighted Articulating Ear Curette™

Bionix® Lighted Articulating Curette combines illumination, magnification and a 90° articulating tip to help healthcare providers take less time to see and safely remove ear canal passage obstructions.



## PROBLEM

Cerumen plugs and foreign bodies in the ear are often difficult to remove. Practitioners lack the ability to visualize or work from behind the occlusion to remove it.

## SOLUTION

The Lighted Articulating Ear Curette brilliantly illuminates the ear canal and provides an articulating tip that allows the practitioner to get behind the occlusion.

## RESULTS

Direct visualization and an articulating ear curette eliminates the blind, top down sweeping motion that can abrade the ear canal improving safety and efficiency.

## PRODUCT OFFERING

**#2511** - Lighted Articulating VersaScoop® (3mm tip) Includes 25 Articulating Ear Curette Tips, 1 Light Source and 1 Magnification Lens.

**#2200** - Extra Light Source

**#2265** - Pack of 5 Magnification Lenses

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

## POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code/Description

**69210** - For potential reimbursement of impacted cerumen.

**69200** - For potential reimbursement of foreign body removal from the ear.

**30300** - For potential reimbursement of foreign body removal from the nose.



EC REP Obelis s.a • Boulevard Général Wahis 53 • 1030 Brussels, BELGIUM  
Tel: +(32) 2.732.59.54 • Fax: +(32) 2.732.60.03 • E-Mail: mail@obelis.net  
Bionix LLC • 1670 Indian Wood Circle, Maumee, Ohio 43537, United States  
Phone: +1 419.727.8421 • Fax: +1 800.455.5678 • www.Bionix.com



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • [www.bionix.com](http://www.bionix.com)

©2021 Bionix LLC | Patent # 8,728,107 | MKT-0033 | Rev. D