

# Lighted Suction

Bionix® Lighted Suction combines illumination and magnification with the efficiency of a suction device, to help healthcare providers take less time to see and safely remove ear canal and nasal passage obstructions.



## PROBLEM

Cerumen plugs, round foreign bodies and soft viscous wax all present problems when it comes to removal. Curettage, and irrigation have often been the solutions but aren't always optimal.

## SOLUTION

Lighted Suction combines the award winning technology of the Lighted Ear Curette™ with Magnification and the simplicity of a suction device.

## RESULTS

A simple and safe suction device that provides complete visualization for suctioning cerumen and round foreign bodies.

## PRODUCT OFFERING

- #2625 - Lighted Suction for Cerumen Removal (1 Light Source, 1 Magnification Lens, 25 Lighted Suction Handles)
- #2630 - Suction Pump
- #2635 - Replacement tubing
- #2200 - Extra Light Source
- #2299 - Light Source Lanyard
- #2265 - Pack of 5 Magnification Lenses

*"I've been using the system (Lighted Suction) for several weeks and have been impressed by the ease of operation and how effective it is compared to curetting cerumen or using an irrigation system."*

- Andrew Schuman, MD

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

## POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code/Description

**69210** - For potential reimbursement of impacted cerumen.

**69200** - For potential reimbursement of foreign body removal from the ear.



**Rx ONLY**

**MADE IN THE U.S.A.**



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • [www.bionix.com](http://www.bionix.com)

©2021 Bionix LLC | US Patent # 9,114,202 | MKT-0037 | Rev. C