



Primary Care



Podiatry



Wound Care



Community Health

CLINIC PLACEMENT



Rapid ABI-Q



Single-Level



TBI w/PPG



Exercise

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



DICOM Compatible

## simpleABI with CUFF-LINK

The simpleABI Automated Cuff-Link System provides tailored testing capabilities and clinical precision for early detection of Peripheral Artery Disease (P.A.D.). It offers versatile testing options, efficiency, and simplifies reimbursement processes, enhancing convenience and safety. With a 2-year warranty, it is the ultimate choice for proactive P.A.D. diagnosis and patient well-being.

### ABI-450CL Exams

- ABI-Q Rapid Screen
- Single-Level
- TBI with PPG
- Exercise

### ABI-450CL Includes:

- Cuffs
  - 2 x 10cm EZ-Clean Vasocuffs
  - 2 x 12cm EZ-Clean Vasocuffs
  - 2 x 2.5cm disposable digit cuffs
- Two-Year Warranty
  - Includes all parts, labor, and accidental breakage



### American-Made Quality

- simpleABI systems are manufactured utilizing the finest quality components and are produced outside of Denver, Colorado.
- Our products are designed to endure years of use, even in the busiest of clinical settings.
- Our industry-leading 2-year warranty provides comprehensive coverage, including labor, parts, and accidental breakage, ensuring complete peace of mind.

### Newman Medical History & Commitment

- Since the 1970s, the Newman family has been at the forefront of ultrasound diagnostics.
- Newman Medical is committed to ongoing product development, guaranteeing the consistent availability of high-quality vascular diagnostic products in the market.
- Newman Medical has a relentless focus on service and establishing long term relationships with our customers

### First-Class Support

- Newman Medical takes great pride in its customer service team, which comprises some of the most seasoned professionals in the industry. On average, our team has over 25 years of experience.
- Our highly proficient field representatives are always ready to offer physicians exceptional support.
- Newman Medical provides physicians with a range of clinical support options, including online and in-person resources, to ensure comprehensive assistance.


## ABI-450CL



### Exams Performed



ABI-Q Rapid Screen



Single-Level



TBI with PPG



Exercise Exams

## Ordering Information

### System

#### SKU

ABI-450CL

#### DESCRIPTION

simpleABI Single-Level Cuff-Link System with Exercise

### Accessories

#### SKU

ABI-120  
STND-500  
ACC-190  
D5  
SPPG  
SFTW-100  
ABI-WARR

#### DESCRIPTION

Microsoft® Windows® PC for use with simpleABI Systems  
Variable-Height Roll Stand for Cuff-Link Systems  
Large simpleABI carry bag  
5Mhz vascular probe only  
PPG Probe for use with Cuff-Link Systems - includes toe cuff  
DICOM software integration  
1 year additional warranty for simpleABI systems

### Replacement Cuffs

#### SKU

CUFF-200  
CUFF-210  
CUFF-210L  
CUFF-220D  
CUFF-240D  
CUFF-SET-242

#### DESCRIPTION

EZ Clean Vascular Cuff, 10cm (ankle, small arm)  
EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)  
EZ Clean Vascular Cuff, 12cm long (above knee)  
Digit Cuff, 2.5cm - Latex-Free  
Digit Cuff, 1.9cm - Latex-Free  
EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



## ABI-450CL Technical Specifications

**Height:** 107 cm (42 inches)

**Weight:** 5.5 kgs (12 pounds)

**Level of Protection against electrical shock:** Type B Applied Part, Class II Equipment

**Comply with the following Standards:** IEC60101-1, IEC60601-2, IEC60601-2-37

**Operating Temperature:** 10° ~ 40°C (50° ~ 104°F)

**Operating Humidity:** 30% ~ 75%

**Transport/Storage Temperature:** -20° ~ 50°C (-4° ~ 122°F)

**Transport/Storage Humidity:** 5% ~ 90%, non-condensing

**Battery Life:** 10 hours

**GTIN:** 00852785008623

## REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

